



puffix
ORTHOPEDICS

CONTENTS

NECK COLLARS

UPPER EXTREMITY

SHOULDER SUPPORTS

LOWER EXTREMITY

PUFFIX KIDS

FOOT HEALTH PRODUCTS

REHABILITATION PRODUCTS

PILLOW AND SEAT CUSHIONS

VARICOSE STOCKINGS

About Us

ATILIM SAĞLIK A.Ş. continues to serve customers with over 400 products and over 50 employees and with superior technology, brand quality and dynamic human resources.

ATILIM SAĞLIK A.Ş. was founded in 2006 in order to establish the standards of a better life by shaping the modern life and pioneering new discoveries in the medical sector.

ATILIM SAĞLIK A.Ş.; which has the best life support brands of its class for its customers as well as its solution partners continues its activities with the aim of becoming a global player in the sector it serves.

ATILIM SAĞLIK A.Ş. It is managed with the aim of contributing to the future of the medical sector to create a progressive and developing future. Behind its success is a management approach centered on customer-focused and efficiency.

As a result of this approach, in 2008, we took our place in reputable trade fair organizations around the world. Our first export was in 2008-2009 to 16 countries; we follow a successful trade policy today by reaching 42 countries.

Located among the pioneers in medical sector in Turkey, ATILIM SAĞLIK A.Ş., which have wide range of services is assessing the best way to reach the target of information and cooperation network.

We continue our pleasant journey with the brand name Puffix, Freeyl, Only Only and Airneb which are representing Turkey in all over the World since 2016



NECK COLLAR



SOFT COLLAR

SUPPORT SOFT COLLAR



Reference Number:

PF - 100

Sizes:

| | |
|---|------------|
| S | 28 - 33 cm |
| M | 34 - 39 cm |
| L | 40 - 45 cm |

Explanation:

It is a soft sponge neck which is designed according to the anatomical neck shape and provides high comfort neck support. Relieves headaches by restricting head and neck movements. The removable outer cover is washable.

Indications:

- Arthritis
- Neck hernia
- Muscle spasms
- Degenerative changes after cervical trauma

Reference Number:

PF - 101

Sizes:

| | |
|---|------------|
| S | 28 - 33 cm |
| M | 34 - 39 cm |
| L | 40 - 45 cm |

Explanation:

It is a soft sponge neck which is designed according to the anatomical neck shape and provides a high comfortable neck support with special support. Relieves headaches by restricting head and neck movements. The removable outer cover is washable.

Indications:

- Arthritis
- Neck hernia
- Muscle spasms
- Degenerative changes after cervical trauma



Reference Number:

PF - 102

Sizes:

| | |
|----|------------|
| S | 32 - 35 cm |
| M | 36 - 39 cm |
| L | 40 - 43 cm |
| XL | 43 - 46 cm |

Explanation:

The anatomically designed 12 mm perforated palastosis is made of foam material. Non-hypoallergenic hypoallergenic foam material does not cause skin irritation. X-ray, CT and MRI are compatible.

Indications:

- Neck hernia
- Soft tissue and cervical muscle strain spasms
- Post-traumatic recovery and protection

Reference Number:

PF - 103

Sizes:

| | |
|----|------------|
| S | 32 - 35 cm |
| M | 36 - 39 cm |
| L | 40 - 43 cm |
| XL | 43 - 46 cm |

Explanation:

The anatomically designed 12 mm perforated palastosis is made of foam material. Non-hypoallergenic hypoallergenic foam material does not cause skin irritation. X-ray, CT and MRI are compatible.

Indications:

- Neck hernia
- Soft tissue and cervical muscle strain spasms
- Post-traumatic recovery and protection

PHILADELPHIA COLLAR



PHILADELPHIA TRACHEOSTOMY COLLAR



Reference Number:

PF - 104

Sizes:

| | |
|---|------------|
| S | 33 - 37 cm |
| M | 38 - 42 cm |
| L | 42 - 46 cm |

Explanation:

Made of 12 mm plastozot foam material. Two plastic supports at the front and rear. Limits flexion, extension, lateral flexion and rotation movements of the neck. Non-hypoallergenic hypoallergenic foam material does not cause skin irritation. X-ray, CT and MRI are compatible.

Indications:

- Stable bone and ligament injuries
- Postoperative protection and stabilization
- Support after Halo use
- Soft tissue injuries
- Suspicious cervical traumas in unconscious patients

Reference Number:

PF - 105

Sizes:

| | |
|---|------------|
| S | 33 - 37 cm |
| M | 38 - 42 cm |
| L | 42 - 46 cm |

Explanation:

Made of 12 mm plastozot foam material. Two plastic supports at the front and rear. Limits flexion, extension, lateral flexion and rotation movements of the neck. Non-hypoallergenic hypoallergenic foam material does not cause skin irritation. X-ray, CT and MRI are compatible.

Indications:

- Stable bone and ligament injuries
- Postoperative protection and stabilization
- Support after Halo use
- Soft tissue injuries
- Suspicious cervical traumas in unconscious patients

SERVICAL COLLAR



SERVICAL COLLAR WITH CHIN



Reference Number:

PF - 106

Sizes:

| | |
|-----|------------|
| S | 32 - 35 cm |
| M | 36 - 39 cm |
| L | 40 - 43 cm |
| XL | 43 - 46 cm |
| XXL | 46 - 50 cm |

Explanation:

It is made of plastic material which is designed according to the anatomical neck shape. Velcro straps can be adjusted to the desired height and provides a long-standing natural posture. Reduces pressure by providing neck stabilization.

Indications:

- Neck hernia
- Whiplash injury
- Post-traumatic protection
- Arthritis

Reference Number:

PF - 107

Sizes:

| | |
|-----|------------|
| S | 32 - 35 cm |
| M | 36 - 39 cm |
| L | 40 - 43 cm |
| XL | 43 - 46 cm |
| XXL | 46 - 50 cm |

Explanation:

It is made of plastic material which is designed according to the anatomical neck shape. Neck support is more restrictive with chin support. Velcro straps can be adjusted to the desired height and provides a long-standing natural posture. Reduces pressure by providing neck stabilization.

Indications:

- Neck hernia
- Arthritis
- Postoperative protection
- Heavy radiculopathy
- Injury on

MINERVA CERVICAL ORTHOSIS



Reference Number:

PF - 108

Sizes:

Neck

| | |
|---|------------|
| S | 33 - 37 cm |
| M | 38 - 42 cm |
| L | 42 - 46 cm |

Shoulder

| | |
|---|--------------|
| S | 70 - 90 cm |
| M | 90 - 110 cm |
| L | 110 - 130 cm |

Explanation:

The outer surface of the polyethylene inner surface is made of terry cloth. It is an orthosis used for upper trocal region stabilization which limits flexion, extension, lateral flexion and rotation movements of the cervical region.

Indications:

- Post-trauma, muscle injuries
- C3 - T2 stable fractures
- Degenetarative joint disorders



ARM SLING (SPONGE)



ARM SLING (FLEXI)



Reference Number:

PF - 200

Sizes:

| | |
|-----|------------|
| S | 30 - 33 cm |
| M | 33 - 36 cm |
| L | 36 - 39 cm |
| XL | 39 - 42 cm |
| XXL | 42 - 45 cm |

Explanation:

The inner surface is made of sponge fabric with a towel. Humerus relaxes the shoulder joint and elbow. Provides shoulder immobilization by fixing the arm to the body. Quality fabric structure provides the comfort of the patient during use.

Indications:

- Post-up immobilization
- Shoulder dislocation and half dislocations
- To prevent the formation of edema in the fracture distal

Reference Number:

PF - 201

Sizes:

| | |
|-----|------------|
| S | 30 - 33 cm |
| M | 33 - 36 cm |
| L | 36 - 39 cm |
| XL | 39 - 42 cm |
| XXL | 42 - 45 cm |

Explanation:

It is made of air perforated fabric. Humerus relaxes the shoulder joint and elbow. Provides shoulder immobilization by fixing the arm to the body. Quality fabric structure provides the comfort of the patient during use.

Indications:

- Post-up immobilization
- Shoulder dislocation and half dislocations
- To prevent the formation of edema in the fracture distal

ACTIVE ARM SLING (30 DEGREE)



Reference Number:

PF - 202

Sizes:

Universal size

Explanation:

The inner surface is made of sponge fabric with a towel. It provides 30 degrees immobilization of the shoulder. Helps mobilization and rehabilitation. Relieves subacromial pressure.

Indications:

- Shoulder endoprosthesis
- Post-up impingement syndrome
- Rotator cuff rupture repair

ACTIVE ARM SLING (45-60 DEGREE)



Reference Number:

PF - 203

Sizes:

Universal size

Explanation:

The inner surface is made of sponge fabric with a towel. It provides immobilization at 45 and 60 degrees abduction in surgical treatment of shoulder, arm and hand. Helps mobilization and rehabilitation. Relieves subacromial pressure.

Indications:

- Shoulder endoprosthesis
- Post-up impingement syndrome
- Rotator cuff rupture repair

SHOULDER SUPPORT

SHOULDER SUPPORT BANDAGE (VELPAU)



Reference Number:

PF - 204

Sizes:

| | |
|-----|------------|
| S | 22 - 24 cm |
| M | 24 - 27 cm |
| L | 27 - 30 cm |
| XL | 30 - 34 cm |
| XXL | 34 - 38 cm |

Explanation:

Made of neoprene material. Increases blood circulation while maintaining body temperature. The compression effect prevents edema, reduces the tone in the muscles, and recovers the pain. The underarm cavity is supported by the pad.

Indications:

- After surgery
- After shoulder prosthesis
- Semi-dislocations
- Therapeutic heat for shoulder and upper arm

Reference Number:

PF - 205

Sizes:

| | |
|-----|------------|
| S | 30 - 33 cm |
| M | 33 - 36 cm |
| L | 36 - 39 cm |
| XL | 39 - 42 cm |
| XXL | 42 - 45 cm |

Explanation:

Made of sponge fabric. It provides immobilization of the shoulder and upper arms in adduction. Bandage wrapped around the shoulder increases stability. Functional design allows the patient to be easily controlled by opening a velcro on the front when needed.

Indications:

- During conservative and post-operative therapy of scapular fractures
- Severe pain in the shoulder joint
- Conservative, post-op proximal humeral head fractures
- Shoulder dislocation and half dislocations
- Shoulder and arm soft tissue injuries

CLAVICLE BANDAGE



ELASTIC ELBOW SUPPORT



Reference Number:

PF - 206

Sizes:

| | |
|----|------------|
| S | 26 - 33 cm |
| M | 34 - 41 cm |
| L | 42 - 49 cm |
| XL | 50 - 57 cm |

Explanation:

The inner surface is made of sponge fabric with a towel. Soft and non-sweating material provides comfort in use with velcro adjustment of the shoulders with the desired tension back and provides smooth posture.

Indications:

- Clavicle and sternum-clavicle dislocation and half dislocations.

Reference Number:

PF - 300

Sizes:

| | |
|-----|------------|
| XS | 20 - 22 cm |
| S | 22 - 24 cm |
| M | 24 - 26 cm |
| L | 26 - 28 cm |
| XL | 28 - 30 cm |
| XXL | 30 - 32 cm |

Explanation:

Made of neoprene material. Elbow helps accelerate the healing process while helping to pass the sensation of pain while maintaining body temperature.

Indications:

- Heat and compression
- Protection against trauma
- Provides protection and support to the elbow
- Protection during and after sports activities.

TENNIS ELBOW SUPPORT WITH EPICONDYLITIS



NEOPREN EPICONDYLITIS BANDAGE



Reference Number:

PF - 302

Sizes:

| | |
|-----|------------|
| XS | 20 - 22 cm |
| S | 22 - 24 cm |
| M | 24 - 26 cm |
| L | 26 - 28 cm |
| XL | 28 - 30 cm |
| XXL | 30 - 32 cm |

Explanation:

Made of neoprene material. There is an anatomically designed silicone pad that provides compression on muscle groups adhering to the lateral epicondyle. Local compression reduces tendon tension and pain.

Indications:

- Lateral epicondylitis (tennis elbow)
- Painful conditions of the elbow
- Increased joint fluid

Reference Number:

PF - 301

Sizes:

Universal size

Explanation:

It is made of neoprene material and it is designed anatomically. A specially designed three-point pad that provides local compression reduces tension in the tendon. The pad is suitable for right and left elbow use.

Indications:

- Medial Epicondylitis (golfer's elbow)
- Lateral Epicondylitis (tennis elbow)

KINETTED ELBOW SUPPORT WITH EPICONDYLITIS

KINETTED ELBOW SUPPORT



Reference Number:

PF - 303

Sizes:

| | |
|-----|------------|
| S | 22 - 24 cm |
| M | 24 - 26 cm |
| L | 26 - 28 cm |
| XL | 28 - 30 cm |
| XXL | 30 - 32 cm |

Explanation:

Manufactured from breathable and flexible mesh fabric. There are silicone pads on both sides of the elbow, which provide compression on the muscle groups that adhere to the condyles. Local compression reduces tendon tension and pain.

Indications:

- Medial epicondylitis
- Lateral epicondylitis
- Post-op or post-traumatic
- Increased joint fluid

Reference Number:

PF - 304

Sizes:

| | |
|-----|------------|
| S | 22 - 24 cm |
| M | 24 - 26 cm |
| L | 26 - 28 cm |
| XL | 28 - 30 cm |
| XXL | 30 - 32 cm |

Explanation:

It is a mesh armrest which is produced according to the anatomic elbow shape which can take the air. It protects body temperature, accelerates blood circulation and healing process. Suitable for use during sports activities.

Indications:

- Post-op or posttraumatic
- Increased joint fluid
- Provides protection and support to the elbow

ELBOW ROM ORTHOSIS



Reference Number:

PF - 305

Sizes:

| | |
|---|------------|
| S | 30 - 36 cm |
| M | 37 - 43 cm |

*There is right and left option.

Explanation:

Anatomical elbow shape designed according to the degree of adjustment of the brace orthosis is easy to adjust the degree adjustment mechanism. The orthosis allows for controlled elbow flexion and extension. The degree of flexion (0-15-30-45-60-75-90-105) and extension (-15-0-15-30-45-60-75-90) can be adjusted. The length of the orthosis can be adjusted according to the arm length with the aluminum support, preventing the fall of the ankle. The shoulder strap provides comfort in use.

Indications:

- In situations where limited and controlled movement is requested on the elbow
- Post-op rehabilitation purposes after elbow endoprosthesis
- Rehabilitation in elbow contracture
- Protective and rehabilitation aid for elbow fractures
- Post-op protection and rehabilitation in elbow ligament instabilities

ELBOW ROM ORTHOSIS



Reference Number:

PF - 306

Sizes:

| | |
|---|------------|
| S | 30 - 36 cm |
| M | 37 - 43 cm |

*There is right and left option.

Explanation:

Anatomical elbow shape designed according to the degree of adjustment of the brace orthosis is easy to adjust the degree adjustment mechanism. The orthosis allows for controlled elbow flexion and extension. The degree of flexion (0-15-30-45-60-75-90-105) and extension (-15-0-15-30-45-60-75-90) can be adjusted.

Indications:

- In situations where limited and controlled movement is requested on the elbow
- Post-op rehabilitation purposes after elbow endoprosthesis
- Rehabilitation in elbow contracture
- Protective and rehabilitation aid for elbow fractures
- Post-op protection and rehabilitation in elbow ligament instabilities

NEOPRENE WRIST SPLINT WITH THUMB SUPPORT



Reference Number:

PF - 400

Sizes:

| | |
|-----|------------|
| S | 14 - 16 cm |
| M | 16 - 18 cm |
| L | 18 - 20 cm |
| XL | 20 - 22 cm |
| XXL | 22 - 24 cm |

Explanation:

Wrist splint made of neoprene material. It provides immobilization of the wrist and thumb in neutral position. It has anatomically shaped aluminum support. It has a bandage that enhances stabilization around the wrist. Neoprene material accelerates blood circulation while maintaining body temperature. The standard is SIZES and it is possible to make personal adjustments with Velcro.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

NEOPRENE WRIST SPLINT



Reference Number:

PF - 401

Sizes:

| | |
|-----|------------|
| S | 14 - 16 cm |
| M | 16 - 18 cm |
| L | 18 - 20 cm |
| XL | 20 - 22 cm |
| XXL | 22 - 24 cm |

Explanation:

Wrist splint made of neoprene material. It provides immobilization of the wrist in neutral position. It has a bandage that enhances stabilization around the wrist. Neoprene material accelerates blood circulation while maintaining body temperature.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

WRIST SPLINT



WRIST SPLINT AIRTEX



Reference Number:

PF - 405

Sizes:

Universal size

*There is right and left option.

Explanation:

It is made of air perforated fabric. It provides immobilization of the wrist in neutral position. It has astigmatically shaped aluminum support. It has a bandage that enhances stabilization around the wrist. Portable aluminum support allows the product to be used for the other wrist.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

Reference Number:

PF - 402

Sizes:

| | |
|---|------------|
| S | 14 - 18 cm |
| M | 18 - 22 cm |

Explanation:

It is made of air perforated fabric. It provides immobilization of the wrist in neutral position. It has a bandage that enhances stabilization around the wrist. Portable aluminum support allows the product to be used for the other wrist

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist

WRIST SPLINT WITH STHUMB SUPPORT AIRTEX



WRIST SPLINT WITH STHUMB SUPPORT



Reference Number:

PF - 403

Sizes:

| | |
|---|------------|
| S | 14 - 18 cm |
| M | 18 - 22 cm |

Explanation:

It is made of air perforated fabric. It provides immobilization of the wrist in neutral position. It has a bandage that enhances stabilization around the wrist. Portable aluminum support allows the product to be used for the other wrist

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist

Reference Number:

PF - 404

Sizes:

Universal size

*There is right and left option.

Explanation:

It is made of air perforated fabric. It provides immobilization of the wrist in neutral position. It has a bandage that enhances stabilization around the wrist. Portable aluminum support allows the product to be used for the other wrist.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

KNITTED WRIST SPLINT



WRIST SPLINT



Reference Number:

PF - 407

Sizes:

| | |
|-----|------------|
| S | 14 - 16 cm |
| M | 16 - 18 cm |
| L | 18 - 20 cm |
| XL | 20 - 22 cm |
| XXL | 22 - 24 cm |

*There is right and left option.

Explanation:

It is made of flexible, non-perspiring, knitted fabric. Supports the immobilization of the wrist in neutral position. The bandage is designed for maximum fit, comfort and effective use, maintains body temperature, accelerates blood circulation and healing process.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

Reference Number:

PF - 406

Sizes:

| | |
|-----|------------|
| S | 14 - 16 cm |
| M | 16 - 18 cm |
| L | 18 - 20 cm |
| XL | 20 - 22 cm |
| XXL | 22 - 24 cm |

*There is right and left option.

Explanation:

It is manufactured from breathable and flexible mesh fabric. It has anatomically shaped aluminum support which allows the wrist to be immobilized in neutral position. Fingers are free from the MCP joint. It has a bandage that provides stabilization around the wrist.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

THERMOPLASTIC RESTING HAND FINGER SPLINT



THERMOPLASTIC ANTI SPASTICITY HAND FINGER SPLINT



Reference Number:

PF - 411

Sizes:

| | |
|---|------------|
| S | 14 - 18 cm |
| M | 18 - 22 cm |

*There is right and left option.

Explanation:

Hand and wrist orthosis produced from thermoplastic material. Holds the wrist and fingers in the rest position. It is used to keep the hand stable in resting position in post-plaster or painful cases.

Indications:

- Plaster after fracture
- In sprains and injuries
- Hemiplegia
- Post-op protection
- Protection from spasticity deformity

Reference Number:

PF - 410

Sizes:

| | |
|---|------------|
| S | 14 - 18 cm |
| M | 18 - 22 cm |

*There is right and left option.

Explanation:

Made of thermoplastic material. It is used in cases of spasticity especially in hand because of the finger design on the splint and resistance to tension.

Indications:

- Spasticity
- Hand - wrist disorders.

THERMOPLASTIC KLEINERT HAND FINGER SPLINT



THERMOPLASTIC DYNAMIC HAND FINGER SPLINT



Reference Number:

PF - 409

Sizes:

| | |
|---|------------|
| S | 14 - 18 cm |
| M | 18 - 22 cm |

*There is right and left option.

Explanation:

Made of thermoplastic material, the inner surface is covered with terry cloth. Allows passive flexion and active extension. It allows early controlled movement after tendon repair in flexor tendon injuries and prevents the formation of adhesions.

Indications:

- Flexor tendon injuries

Reference Number:

PF - 408

Sizes:

| | |
|---|------------|
| S | 14 - 18 cm |
| M | 18 - 22 cm |

*There is right and left option.


Explanation:

It is made of thermoplastic material and its inner surface is covered with towel fabric. Passively extends wrist and fingers. Allows active finger flexion. It allows early controlled movement after tendon repair in extensor tendon injuries and prevents the formation of adhesions.

Indications:

- Extensor tendon injuries

BASEBALL SPLINT - FROG SPLINT - FINGER STABILIZER SPLINT - STRAIGHT FINGER SPLINT



Reference Number:

PF - 502
 PF - 504
 PF - 503
 PF - 500
 PF - 501

Sizes:

| | |
|---|------------|
| S | 4,5 - 5 cm |
| M | 5 - 6 cm |
| L | 6 - 7 cm |

Explanation:

Special aluminum, the inner surface is covered with plastozot finger splint. It provides the immobilization of the finger in extension.

Indications:

- Finger fractures
- Tendon and nerve injuries
- Finger injuries and sprains

Reference Number:

PF - 501

Sizes:

| | |
|-----|------------|
| 1 | 4 - 4.5 cm |
| 2 | 4.5 - 5 cm |
| 3 | 5 - 5.5 cm |
| 4 | 5.5 - 6 cm |
| 5 | 6 - 6.5 cm |
| 5.5 | 7.5 - 8 cm |
| 6 | 8 - 8.5 cm |
| 7 | 8.5 - 9 cm |

Explanation:

It is produced from special plastic material. The finger is free from the PIP joint.

Indications:

- Mallet finger deformity
- PIP joint fractures, sprains and injuries



CHEST CORSET WITH STEEL SUPPORT

DORSOLUMBAR CORSET



Reference Number:

PF - 600

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. The steel supports in the sternal region reduce the pain of movement and support the sternum by providing compression to the rib cage, giving the patient a feeling of confidence and comfort. It helps to provide protection against pain caused by sudden movement.

Indications:

- Fractures of costa and sternum
- Thoracic surgery and post-traumatic
- Support post-surgical sternum

Reference Number:

PF - 601

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. The shoulder straps that provide postural support are crossed under the armpit and crossed at the rear and joined with velcro straps. The crossbands provide extra support in the lumbar region. In the back, 2 steel bales that can be attached and removed from both sides of the spine are inclined according to the physiological curvatures of the spine and provide excellent stability in the body. It gives a feeling of confidence and comfort to the patient in daily life.

Indications:

- Osteoporosis
- Pre and postoperative support
- Vertebral traumas
- Mild kyphosis

ABDOMINAL CORSET



Reference Number:

PF - 602

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. It gives support to the abdomen with its flexible structure. Can be adjusted with the desired tightness of velcro. Reduces the sensation of pain with movement in daily life. It can be worn easily under clothing.

Indications:

- Abdominal surgery and postnatal
- Abdominal sagging
- After umbilical and ventral hernia operations

BELLY FINGER CORSET



Reference Number:

PF - 604

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

Made of flexible fabric. There are compression pads applied to the umbilical hernia. Designed to push back the herniated tissue does not interfere with the bandage daily. It can be worn easily under clothing.

Indications:

- Umbilical hernia

LUMBOSTAD CORSET (26CM) WITH STRAP



LUMBOSACRAL CORSET (32CM) WITH STRAP



Reference Number:

PF - 603

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. It is 26 cm in height and supports up to the lower back and lower back. There are two plastic bumps in the front and rear and 2 steel bumps on the back of the spine. The cross band belt provides extra stabilization in the lumbar region. This reduces pain by reducing the load on the discs.

Indications:

- Immobilization before and after surgical procedure
- Lumbago, degenerative disc lesions
- Acute or chronic low back pain
- Spondylolysis, spondylolisthesis

Reference Number:

PF - 605

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. Its height is 32 cm and supports up to the waist and lower back. There are two plastic bumps in the front and rear and 2 steel bumps on the back of the spine. The cross band belt provides extra stabilization in the lumbar region. This reduces pain by reducing the load on the discs.

Indications:

- Immobilization before and after surgical procedure
- Lumbago, degenerative disc lesions
- Acute or chronic low back pain
- Spondylolysis, spondylolisthesis

POSTUREX CORSET

KNITTED LUMBOSACRAL CORSET



Reference Number:

PF - 606

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. On both sides of the spine there are plastic blades that can be removed and removed. The flexible shoulder straps are crossed under the armpit and crossed to the rear and fastened with velcro straps. Cross bands provide sensory feedback that prevents bad posture disorder by providing postural support. It can be worn easily under clothing.

Indications:

- Posture disorders
- Back and shoulder muscles weakness

Reference Number:

PF - 607

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

High air permeability, which prevents sweat is specially designed corset back with the 4 steel mesh ballet to the anatomical body shape. Steel bales can be manually tilted according to lumbar lordosis. Its conformity to the anatomical structure and fabric structure provides a perfect fit to the body. Lumbar area relieves the tension in the muscles, relieves pain, gives a feeling of comfort.

Indications:

- Post-op protection and immobilization
- Lumbago, disc lesions
- Spondylolysis, spondylolisthesis

KNITTED LUMBOSACRAL CORSET WITD PAD



WOOLEN CORSET CAMEL



Reference Number:

PF - 608

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

High air permeability, which prevents sweat is specially designed corset back with the 4 steel mesh ballet to the anatomical body shape. Steel bales can be manually tilted according to lumbar lordosis. Its conformity to the anatomical structure and fabric structure provides a perfect fit to the body. Lumbar area relieves the tension in the muscles, relieves pain, gives a feeling of comfort.

Indications:

- Post-op protection and immobilization
- Lumbago, disc lesions
- Spondylolysis, spondylolisthesis

Reference Number:

PF - 609

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

Made of wool fabric. Provides temperature and compression while maintaining body temperature in waist area.

Indications:

- Muscle spasms
- Rheumatic and sciatic pain
- To keep the abdomen, waist and kidneys warm

PREGNANCY CORSET



Reference Number:

PF - 610

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. The abdomen can be used before and after prenatal support. Pregnancy line supports belly and prevents back and back pain. It gives the person a sense of confidence and comfort in daily life movements. After the delivery, the large part of the corset is brought to the front and is used to collect the sagging of the abdomen.

Indications:

- Back and low back pain during pregnancy
- Postpartum abdomen recovery

COLOSTOMY CORSET



Reference Number:

PF - 611

Sizes:

| | |
|-----|--------------|
| S | 75 - 85 cm |
| M | 85 - 95 cm |
| L | 95 - 105 cm |
| XL | 105 - 115 cm |
| XXL | 115 - 125 cm |

Explanation:

It is made of flexible and non-sweaty fabric. With the 2 plastic ballet on both sides of the flexible fabric and stomatal opening, preventing the possible leakage that may occur during daily activities in the stoma bag. It gives the patient a sense of trust and peace. It is easily used under clothing.

Indications:

- Protect the stoma bag
- Support after ostomy surgery
- Reduce the risk of parastomal hernia after ostomy surgery.

T.L.S.O. CORSET

L.S.O. CORSET



Reference Number:

PF - 701

Sizes:

| | |
|---|--------------|
| S | 60 - 85 cm |
| M | 85 - 110 cm |
| L | 110 - 140 cm |

Explanation:

Inside the body brace which is made of sponge cloth towels. The corset is provided with a 3 mm support on the abdomen and waist, providing a rigid support and 2 mm on the sternum. By limiting spine movements to a certain extent, it improves anatomically correct posture and improves stability.

Indications:

- After surgery
- Postop fusion
- Compression fractures
- Fractures in the spine

Reference Number:

PF - 700

Sizes:

| | |
|---|--------------|
| S | 60 - 85 cm |
| M | 85 - 110 cm |
| L | 110 - 140 cm |

Explanation:

The interior of the towel is made of sponge cloth lumbo-sacral orthosis. The corset has anatomically designed, 3 mm plastic supports that provide a rigid support on the abdomen and waist. By limiting the lumbar vertebra movements to a certain extent, it provides anatomically correct posture and improves stability.

Indications:

- Stable compression fractures
- Post-traumatic lumbar stabilization
- Postoperative immobilization of lumbar vertebra fractures
- Lumbo-sacral facet syndrome

PECTUS ORTHOSIS



HYPEREXSTANCE CORSET



Reference Number:

PF - 704

Sizes:

Universal size

Explanation:

The inner part of the sponge is made of sponge fabric is an orthosis. The corset is provided with a 3 mm support in the abdomen and waist, providing a rigid support and 2 mm in the sternum. By limiting spine movements to a certain extent, it improves anatomically correct posture and improves stability.

Indications:

- After surgery
- Post-op fusion
- Compression fractures
- Fractures in the spine

Reference Number:

PF - 703

Sizes:

| | |
|---|------------|
| S | 40 - 49 cm |
| M | 50 - 60 cm |

Explanation:

Made of special aluminum, skin contact surfaces covered with a towel pad, 3 point principle of the body orthosis is appropriate. It is a body orthosis that restricts the flexion of the trunk and promotes the position of hyperextension. With the brace length and width adjustment may be threaded system. By limiting the movements of the spine to a certain extent, it provides an anatomically correct posture and improves the stability and accelerates the healing process.

Indications:

- Osteoporosis compression fractures
- Kyphosis
- Compression fractures

HYPEREXSTANCE CORSET



Reference Number:

PF - 702

Sizes:

| | |
|---|------------|
| S | 40 - 46 cm |
| M | 47 - 53 cm |
| L | 54 - 60 cm |

Explanation:

Made of special aluminum, the skin-contacting surfaces are leather-covered body orthosis. Length and width adjustment can be made on the orthotics by means of the screw system on it. The pelvic bar is a movable system that allows you to sit and relax. Sternal pads reduce the regional pressure and provide comfortable operation.

Indications:

- Post-operative orthotic care and bracing
- Post-op instable vertebral body fractures
- Stable compression fractures of the lumbar and lower thoracic vertebrae
- After intervertebral disc surgery

INGUINAL TRUSS (SINGLE SIDED)

INGUINAL TRUSS (DOUBLE SIDED)



Reference Number:

PF - 001

Sizes:

| | |
|---|--------------|
| S | 75 - 90 cm |
| M | 90 - 105 cm |
| L | 105 - 120 cm |

*There is right and left option.

Explanation:

Made of flexible cotton fabric. Specially designed anatomic pad and bandage to push back the herniated tissue does not interfere with normal daily activities. It is used to prevent recurrences after surgery and to support healing. It provides the perfect solution for support and comfort to prevent pain and discomfort. The buckled system provides comfort in use.

Indications:

- Single-sided inguinal hernia

Reference Number:

PF - 002

Sizes:

| | |
|---|--------------|
| S | 75 - 90 cm |
| M | 90 - 105 cm |
| L | 105 - 120 cm |

Explanation:

Made of flexible cotton fabric. Specially designed anatomic pad and bandage to push back the herniated tissue does not interfere with normal daily activities. It is used to prevent recurrences after surgery and to support healing. It provides the perfect solution for support and comfort to prevent pain and discomfort. The buckled system provides comfort in use.

Indications:

- Double-sided inguinal hernia



SIMPLE KNEE SUPPORT



PATELLA KNEE SUPPORT



Reference Number:

PF - 800

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

Made of neoprene material. Supports the knee joint, alleviates the burden on the insertion and reduces the pain. Neoprene material, while maintaining the body temperature increases blood circulation, the compression effect prevents edema, reduces the tone in the muscles, relieves pain.

Indications:

- Post-operative protection and rehabilitation
- Protection during sports activities
- Increased joint fluid
- Osteoarthritis

Reference Number:

PF - 801

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. The circumference of the patella is open and supported by the pad. It is ideal for patients with pain in the patella and having problems after surgery.

Indications:

- Post-operative protection and rehabilitation
- Protection during sports activities
- Increased joint fluid
- Osteoarthritis

PATELLA AND LIGAMENT KNEE SUPPORT



Reference Number:

PF - 802

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. The circumference of the patella is open and supported by the pad. The flexible brackets on the side increase the stabilization and thus reduce the string load. Neoprene material accelerates blood circulation while maintaining body temperature.

Indications:

- Meniscus problems
- Arthritis, bursitis
- Pain from knee cap
- Postoperative protection
- 1st degree sprain and ligament instabilities

HINGED KNEE SUPPORT



Reference Number:

PF - 803

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. Aluminum joints on both sides allow maximum controlled support for directional movement. Prevents index hyperextension. It is recommended during rehabilitation and treatment after severe surgeries and knee injuries.

Indications:

- Postoperative protection and immobilization
- Meniscus lesions
- Ligament injuries and instability
- 1st and 2nd degree sprains and strains

KNEE SUPPORT FOR ANT.

KNEE SUPPORT FOR ANT. CRUCIATE LIGAMENT



Reference Number:

PF - 804

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. The patella area is open and supported by pad. The medial and lateral spiral balines increase the stability around the joint. It provides stabilization of the anterior cruciate ligaments and patella with the crossbands where the three-point principle is applied.

Indications:

- Meniscus lesions
- Postoperative protection and stabilization
- Anterior cruciate ligament injuries and knee instability
- Patellar stabilization

Reference Number:

PF - 805

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. The circumference of the patella is open and supported by the pad. The medial and lateral spiral balines increase the stability around the joint. It provides stabilization of the anterior cruciate ligaments and patella with the crossbands where the three-point principle is applied.

Indications:

- Meniscus lesions
- Postoperative protection and stabilization
- Anterior cruciate ligament injuries and knee instabilities
- Patellar stabilization

PATELLAR TENDON BAND



Reference Number:

PF - 806

Sizes:

Universal size

Explanation:

It is made of neoprene material and it is designed anatomically. There is an anatomically designed pad that supports the patellar tendon. The compression effect reduces tension and pain in the tendon. Neoprene is flexible and comfortable to use with material content. Both strings can be applied.

Indications:

- Patellar tendinitis
- Patella femoral pain
- Chondromalacia patella
- Os-good Schlatter's disease

WOOLEN KNEE SUPPORT CAMEL



Reference Number:

PF - 807

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

Made of wool fabric. It provides temperature and compression while maintaining the body temperature around the patella.

Indications:

- Arthritis
- Temperature purpose

KNITTED SIMPLE KNEE SUPPORT



Reference Number:

PF - 808

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of breathable, flexible, non-sweat-proof mesh fabric. The heat effect it provides in the knee joint provides long-term relief on the aching joints, muscles and surrounding tissues. Since it does not press the knee cap, it can be preferred in daily use.

Indications:

- Osteoarthritis
- Increased joint fluid
- Sprains
- Protection during sports activities

KNEE PATELLA KNEE SUPPORT



Reference Number:

PF - 809

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of flexible, non-perspiring, knitted fabric. High quality woven fabric provides comfortable use throughout the day. The silicone pad surrounding the patella provides support around the knee cap. The medial and lateral spiral balines increase the stability around the joint. Supports the knee joint, alleviates the burden on the insertion and reduces the pain.

Indications:

- Osteoarthritis, bursitis
- Increased joint fluid
- Knee cap-induced pain
- Postoperative protection and stabilization

KNEE PATELLA AND LIGAMENT KNEE SUPPORT



Reference Number:

PF - 810

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of flexible, non-per-spiring, knitted fabric. High quality woven fabric provides comfortable use throughout the day. The silicone pad surrounding the patella provides support around the knee cap. The medial and lateral spiral balines increase the stability around the joint. Supports the knee joint, alleviates the burden on the insertion and reduces the pain.

Indications:

- Meniscus problems
- Arthritis, bursitis
- Pain from knee cap
- Postoperative protection
- 1st degree sprain and ligament instabilities

KNITTED HINGED KNEE SUPPORT



Reference Number:

PF - 811

Sizes:

| | |
|-----|------------|
| S | 31 - 34 cm |
| M | 34 - 37 cm |
| L | 37 - 40 cm |
| XL | 40 - 44 cm |
| XXL | 44 - 48 cm |

Explanation:

It is made of flexible, non-per-spiring, knitted fabric. Aluminum joints on both sides allow maximum controlled support for directional movement. The silicone pad around the patella reduces pressure by protecting the knee cap and prevents index hyperextension. It is recommended during rehabilitation and treatment after severe surgeries and knee injuries.

Indications:

- Postoperative protection and immobilization
- Meniscus tears
- Ligament injuries and instability
- 1st and 2nd degree buckling and stresses
- Ligament Injuries



Reference Number:

PF - 812

Sizes:

| | |
|---|-------|
| 1 | 45 cm |
| 2 | 55 cm |
| 3 | 65 cm |

Explanation:

The inner surface of which is made of sponge cloth towels. Provides immobilization of the index extension. Removable aluminum supports are available that provide posterior, medial and lateral stabilization. It is supportive for knee injuries and postoperative surgery. Standard SIZES are available in lengths of 45 cm, 55 cm and 65 cm. With a three-piece pad design, the legs can be easily adjusted to the environment.

Indications:

- Meniscal injury after surgery
- Acute knee sprains and injuries
- Helps post-operative rest, protection and immobilization
- Protection after ligament injuries
- Osteochondral fractures
- Knee joint traumas

Reference Number:

PF - 813

Sizes:

| | |
|-------|-------|
| SHORT | 47 cm |
| LONG | 58 cm |

Explanation:

Anatomically designed aluminum bar and skin contact surfaces are terry cloth, which is a highly adjustable knee orthosis. Flexion (0 °, 15 °, 30 °, 45 °, 60 °, 75 °, 90 °, 105 °, 120 °) and extension (-15 °, 0 °, 15 °, 30 °, 45 °, 60 °, 75 °, 90 °) degree adjustment. It offers excellent stabilization and support after knee injuries. Surgery and post-trauma rehabilitation period is used as a preservative and R. O. directory provides control, 2 lengths.

Indications:

- Post-operative grade mobilization
- After surgery on the meniscus injury
- 4. High tibial osteotomy
- Ambulation (mobilization)
- Post-operative controlled and graded joint mobilization in ligament injuries

ANGLE - ADJUSTED BRACE (PUFFIX PLUS)

HIP BRACE



Reference Number:

PF - 814

Sizes:

Universal size

Explanation:

Anatomically designed aluminum bar and skin contact surfaces are terry cloth, which is a highly adjustable knee orthosis. Flexion (0°, 15°, 30°, 45°, 60°, 75°, 90°, 105°, 120°) and extension (-15°, 0°, 15°, 30°, 45°, 60°, 75°, 90°) degree adjustment. It offers excellent stabilization and support after knee injuries. It is used as a preservative during surgery and post traumatic rehabilitation period and provides control of the index R.O.M, it is one size. Size adjustment can be made with the buttons on the product.

Indications:

- Post-operative grade mobilization
- After surgery on the meniscus injury
- 4. High tibial osteotomy
- Ambulation (mobilization)
- Post-operative controlled and graded joint mobilization in ligament injuries
- Stable femoral fractures

Reference Number:

PF - 003

Sizes:

| Calf | Waist Size |
|-----------------|--------------|
| S : 40 - 45 cm | 60 - 75 cm |
| M : 45 - 52 cm | 75 - 90 cm |
| L : 52 - 60 cm | 90 - 105 cm |
| XL : 60 - 70 cm | 105 - 120 cm |

Explanation:

Anatomically designed is made of thermoplastic material. Inner surface that contacts the skin is covered with cloth towel. It is the hip abduction orthosis that provides stabilization of the hip joint. Up to 120 degrees of flexion can be adjusted in 15-degree increments. Abduction angle can be adjusted.

Indications:

- Stabilization after total hip replacement or revision
- Femoral head resection
- Repeated cases of hip dysplasia
- After femoral or acetabular osteotomy

SIMPLE ANKLE SUPPORT



Reference Number:

PF - 555

Sizes:

| | |
|-----|------------|
| S | 21 - 23 cm |
| M | 23 - 25 cm |
| L | 25 - 27 cm |
| XL | 27 - 29 cm |
| XXL | 29 - 31 cm |

Explanation:

It is made of neoprene material, and the anatomically designed. Supports the ankle, reducing pain. Neoprene material accelerates blood circulation while maintaining body temperature.

Indications:

- Protection against minor fractures
- 1st degree sprains and strains
- Post-buckling stabilization and protection
- Compression edema after injury

Reference Number:

PF - 556

Sizes:

| | |
|-----|------------|
| S | 21 - 23 cm |
| M | 23 - 25 cm |
| L | 25 - 27 cm |
| XL | 27 - 29 cm |
| XXL | 29 - 31 cm |

Explanation:

It manufactured from Neoprene material. There are pads on both sides of the ankle supporting internal and external malleolus. The eight-fold elastic band provides extra protection. It supports the ankle, relieves pain, reduces edema.

Indications:

- After surgery and injuries protective and support for
- After Gypsum
- In cases of edema, hematoma and effusion
- Mild ligament injuries and instability
- Joint inflammations and calcifications

MALLEOL ANKLE SUPPORT



LIGAMENT ANKLE SUPPORT



Reference Number:

PF - 557

Sizes:

| | |
|-----|------------|
| S | 21 - 23 cm |
| M | 23 - 25 cm |
| L | 25 - 27 cm |
| XL | 27 - 29 cm |
| XXL | 29 - 31 cm |

Explanation:

It is manufactured from Neoprene material. The inner and outer plastic plates, which provide stabilization around the malleoli, are manufactured separately as internal and external. It can be opened completely and provides comfort in use. The elastic cross band that wraps in eight forms is an extra protection bandage. Apply the compression effect that the ankle support, reduces edema and pain.

Indications:

- Protection or support after surgery and injury
- Bond instability and tears
- Protective and edema preventive

KNITTED SIMPLE ANKLE SUPPORT



Reference Number:

PF - 558

Sizes:

| | |
|-----|------------|
| S | 19 - 21 cm |
| M | 21 - 23 cm |
| L | 23 - 25 cm |
| XL | 25 - 27 cm |
| XXL | 27 - 29 cm |

Explanation:

It is made of flexible, non-perspiring, knitted fabric. The effect of heat on the ankle provides long-term relief on the aching joints, muscles and surrounding tissues. The compression effect reduces edema and pain.

Indications:

- Postoperative protection
- After minor fractures
- Protection against sprains
- 1st degree sprains and strains

KNITTED MALLEOL ANKLE SUPPORT



KNITTED ACHILLES TENDON ANKLE SUPPORT



Reference Number:

PF - 559

Sizes:

| | |
|-----|------------|
| S | 19 - 21 cm |
| M | 21 - 23 cm |
| L | 23 - 25 cm |
| XL | 25 - 27 cm |
| XXL | 27 - 29 cm |

Explanation:

It is made of flexible, non-perspiring, knitted fabric. There are anatomically silicone pads around the malleoles. Supporting the ankle, the compression effect reduces pain and edema.

Indications:

- Protective and supportive purpose after surgery and injury
- In cases of edema, hematoma and effusion
- Severe ligament injuries and instability
- Joint inflammation and calcification
- Early mobilization aid

Reference Number:

PF - 560

Sizes:

| | |
|-----|------------|
| S | 19 - 21 cm |
| M | 21 - 23 cm |
| L | 23 - 25 cm |
| XL | 25 - 27 cm |
| XXL | 27 - 29 cm |

Explanation:

Alt is made of flexible, non-perspiring, knitted fabric. Silicone support is available on both sides of the Achilles tendon and under the heel. There are also portable heels placed under the other heel to prevent imbalance between both legs. It protects the heel and the Achilles tendon and reduces the tension in the tendon.

Indications:

- Acute and chronic inflammation of the Achilles tendon
- Achilles tendon pain
- After surgery or injury

ANKLE SUPPORT WITH PAD



ANKLE SUPPORT WITH AIR



Reference Number:

PF - 561

Sizes:

Universal size

Explanation:

The padded ankle is the last foot ankle with a sponge pad designed for the inner and outer stabilization of the ankle without hindering the natural joint movement. It is ideal for the treatment and rehabilitation of mild to moderate sprains of the ankle.

Indications:

- Post-plaster stabilization
- Post-sprain compression and support
- Pre-op and Post-op ligament injuries
- Chronic instability of the ankle

Reference Number:

PF - 562

Sizes:

Universal size

Explanation:

Pneumatic ankle is a foot pad with air pad designed for the shape of the inner and outer stabilization of the ankle without hindering the movement of the natural joint. It is ideal for the treatment and rehabilitation of mild to moderate sprains of the ankle. Thanks to the inflatable air pad, stabilization adjustment can be made at the desired level.

Indications:

- Post-plaster stabilization
- Post-sprain compression and support
- Pre-op and Post-op ligament injuries
- Chronic instability of the ankle

DORSI - FLEXION BANDAGE



REFLEKS AFO



Reference Number:

PF - 563

Sizes:

Universal size

Explanation:

It is manufactured from Neoprene material. Foot ankle used in low foot deformity which keeps the ankle dorsiflexed at desired tension. Provides a visible improvement in the patient's walk. Available in shoes. Standard size and can be used for right and left feet.

Indications:

- Muscle weakness in paralyzed patients
- Low foot deformity
- Weak dorsiflexion

Reference Number:

PF - 564

Sizes:

| | |
|---|----------------|
| S | 34 - 37 Number |
| M | 38 - 41 Number |
| L | 42 - 45 Number |

Explanation:

It is an ankle orthosis made of plastic material. Increases lateral stability by providing static dorsiflexion of the wrist. Velcro strap around the ankle provides comfort in use. It can be shaped with heat gun when necessary.

Indications:

- In cases where the ankle is required to be kept at 90 degrees
- Low foot deformity that develops when nerve injuries such as hemiplegia and lumbar hernia occur.

FOOT GUARD ADJUSTABLE



Reference Number:

PF - 566

Sizes:

| | |
|---|---------|
| S | 35 - 37 |
| M | 38 - 40 |
| L | 41 - 46 |

Explanation:

Thermo plastic material terry fabric is covered and removed.
T-band ankle orthosis. 90 degrees of the ankle and T-band prevents rotation. Foot and ankle in multi position Provides rest. Eskar wounds It prevents. angles securing T Prevents rotation with tape. It is preferred in hospitalized patients.

Indications:

- To maintain the position of the foot in the inpatient

FOOT GUARD



Reference Number:

PF - 565

Sizes:

| | |
|---|----------------|
| M | 35 - 39 Number |
| L | 40 - 44 Number |

Explanation:

It is made of thermoplastic material. Prevents hip rotation by fixing the ankle at 90 degrees.

Indications:

- To maintain the position of the foot in the inpatient

R.O.M. WALKER



ACHILLES BOOT



Reference Number:

PF - 567

Sizes:

| | |
|---|----------------|
| S | 35 - 39 Numara |
| M | 40 - 43 Numara |
| L | 44 - 47 Numara |

Explanation:

Ankle dorsi and plantar flexion degree can be adjusted to help early mobilization of ankle orthosis. The plantar and dorsi flexion of the ankle can be adjusted to 0 ° -7.5 ° -15 ° -22.5 ° -30 ° -37.5 ° -45 degrees. The rocker model design at the base supports the natural walking pattern, allowing you to walk more comfortably.

Indications:

- Distal tibia and fibula fractures in stable distal tibia and fibula fractures or internal fixator
- Achilles tendon ruptures
- 5. Stable foot and ankle fractures
- Postoperative stabilization for tendon, ligament and soft tissue injuries

Reference Number:

PF - 568

Sizes:

| | |
|---|----------------|
| S | 35 - 39 Numara |
| M | 40 - 43 Numara |
| L | 44 - 47 Numara |

Explanation:

The inner surface is made of sponge fabric. An ankle orthosis that provides immobilization of the ankle. The rocker model design at the base part supports the natural walking pattern and provides more comfortable walking.

Indications:

- Distal tibia and fibula fractures in stable distal tibia and fibula fractures or internal fixator
- Achilles tendon ruptures
- 5. Stable foot and ankle fractures
- Postoperative stabilization for tendon, ligament and soft tissue injuries

STANDARD SHOULDER SUPPORT



Reference Number:

PF - 900

Sizes:

Universal size

Explanation:

The inner surface is made of sponge fabric with a towel. Increases blood circulation while maintaining body temperature. The compression effect prevents edema, reduces the tone in the muscles, and recovers the pain. The underarm cavity is supported by the pad. Standard size and can be personalized with velcro.

Indications:

- After surgery
- After shoulder prosthesis
- Semi-dislocations
- Therapeutic heat for shoulder and upper arm

STANDARD VELPAU BANDAGE



Reference Number:

PF - 901

Sizes:

Universal size

Explanation:

Made of sponge fabric. It provides immobilization of the shoulder and upper arms in adduction. Bandage wrapped around the shoulder increases stability. Functional design allows the patient to be controlled easily by opening a velcro on the front when needed.

Indications:

- During conservative and post-operative therapy of scapular fractures
- Severe pain in the shoulder joint
- Conservative, post-op proximal humeral head fractures
- Shoulder dislocation and half dislocations
- Shoulder and arm soft tissue injuries

STANDARD CLAVICLE BANDAGE



Reference Number:

PF - 902

Sizes:

Universal size

Explanation:

The inner surface is made of sponge fabric with a towel. Soft and non-sweating material provides comfort in use with velcro adjustment of the shoulders with the desired tension back and provides smooth posture. Standard size and can be personalized with velcro.

Indications:

- Clavicle and sternum-clavicle dislocation and half dislocations.
- Clavicle fractures

STANDARD ELBOW SUPPORT



Reference Number:

PF - 903

Sizes:

Universal size

Explanation:

Made of neoprene material. Elbow helps accelerate the healing process while helping to pass the sensation of pain while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- Heat and compression
- Protection against trauma
- Provides protection and support to the elbow
- Protection during and after sports activities.

STANDARD ELBOW SUPPORT WITH EPICONDYLITIS



Reference Number:

PF - 904

Sizes:

Universal size

Explanation:

Made of neoprene material. Bandages in the form of wrapping the hands and provides protection. Helps protect the heat from wrist injuries, sprains and traumas. Bandage can also help prevent injury and dislocations during sports and physical activities.

Indications:

- Wrist injury and sprains
- Protection during sports activities

STANDARD DE QUERVAIN SPLINT



Reference Number:

PF - 907

Sizes:

Universal size

Explanation:

Neoprene material is a thumb fingerboard. It has aluminum support that can be adjusted with a slope that restricts thumb movements. Velcro straps provide a comfortable pressure. Neoprene material accelerates blood circulation while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- De Quervain syndrome
- Injury of the thumb ligament
- 1.CMC and MCP joint irritation

STANDARD WRIST SPLINT



Reference Number:

PF - 906

Sizes:

Universal size

Explanation:

Wrist splint made of neoprene material. It provides immobilization of the wrist in neutral position. It has a bandage that enhances stabilization around the wrist. Neoprene material accelerates blood circulation while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful conditions of the wrist
- Plaster after fracture
- Injury and after surgery

STANDARD WRIST SPLINT WITH THUMB SUPPORT



Reference Number:

PF - 905

Sizes:

Universal size

Explanation:

Wrist splint made of neoprene material. It provides immobilization of the wrist and thumb in neutral position. It has an anatomically shaped aluminum support. It has a bandage that enhances stabilization around the wrist. Neoprene material accelerates blood circulation while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

STANDARD ABDOMINAL CORSET



Reference Number:

PF - 908

Sizes:

Universal size

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. It gives support to the abdomen with its flexible structure. Can be adjusted with the desired tightness of velcro. Reduces the sensation of pain with movement in daily life. It can be worn easily under clothing. Standard size and can be personalized with velcro.

Indications:

- Abdominal surgery and postnatal
- Abdominal sagging
- After umbilical and ventral hernia operations

STANDARD LUMBOSTAD CORSET (32CM)



Reference Number:

PF - 909

Sizes:

Universal size

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. The height is 32cm, up to the waist and bottom backpack. There are two plastic bumps in the front and rear and two steel bumps that can be removed and removed at the back of the keel. Cross-belt belt provides extra stabilization of lumbar setup. In this way, by reducing the load on the discs to clear the pain. Standard size and personal adjustment with velcro.

Indications:

- Immobilization before and after surgical procedure
- Lumbago, degenerative disc lesions
- Acute or chronic low back pain
- Spondylolysis, spondylolisthesis

STANDARD LUMBOSTAD CORSET (26CM)



Reference Number:

PF - 910

Sizes:

Universal size

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. It is 26 cm in height and supports up to the lower back and lower back. There are two plastic bumps in the front and rear and 2 steel bumps on the back of the spine. The cross band belt provides extra stabilization in the lumbar region. This reduces pain by reducing the load on the discs. Standard size and can be personalized with velcro.

Indications:

- Immobilization before and after surgical procedure
- Lumbago, degenerative disc lesions
- Acute or chronic low back pain
- Spondylolysis, spondylolisthesis

STANDARD BASIC KNEE SUPPORT



Reference Number:

PF - 911

Sizes:

Universal size

Explanation:

Wrist splint made of neoprene material. It provides immobilization of the wrist and thumb in neutral position. It has anatomically shaped aluminum support. It has a bandage that enhances stabilization around the wrist. Neoprene material accelerates blood circulation while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- Post-operative protection and rehabilitation
- Protection during sports activities
- Increased joint fluid
- Osteoarthritis

STANDARD PATELLA KNEE SUPPORT



Reference Number:

PF - 912

Sizes:

Universal size

Explanation:

It is made of neoprene material and it is designed anatomically. The circumference of the patella is open and supported by the pad. It is ideal for patients with pain in the patella and having problems after surgery. Standard size and can be personalized with velcro.

Indications:

- Post-operative protection and rehabilitation
- Protection during sports activities
- Increased joint fluid
- Osteoarthritis

STANDARD PATELLA AND LIGAMENT KNEE SUPPORT



Reference Number:

PF - 913

Sizes:

Universal size

Explanation:

It is made of neoprene material and it is designed anatomically. The circumference of the patella is open and supported by the pad. The flexible brackets on the side increase the stabilization and thus reduce the string load. Neoprene material accelerates blood circulation while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- Meniscus problems
- Arthritis, bursitis
- Pain from knee cap
- Postoperative protection
- 1st degree sprain and ligament instabilities

STANDARD HINGED KNEE SUPPORT



Reference Number:

PF - 914

Sizes:

Universal size

Explanation:

It is made of neoprene material and it is designed anatomically. Aluminum joints on both sides allow maximum controlled support for directional movement. Prevents index hyperextension. It is recommended during rehabilitation and treatment after severe surgeries and knee injuries. Standard size and can be personalized with velcro.

Indications:

- Postoperative protection and immobilization
- Meniscus lesions
- Ligament injuries and instability
- 1st and 2nd degree sprains and strains

STANDARD THIGH SUPPORT



Reference Number:

PF - 916

Sizes:

Universal size

Explanation:

Made of neoprene material. Pain by reducing muscle tone in the upper thigh supports preserving the hafifletir. baldırlık body heat while helping to pass the sensation of pain, accelerates the healing process. Standard size and can be personalized with velcro.

Indications:

- Upper calf muscle stresses protection during sports activities
- The upper calf is painful
- The komperasyo

STANDARD CALF SUPPORT



Reference Number:

PF - 915

Sizes:

Universal size

Explanation:

Made of neoprene material. Supports the lower calf by reducing the tone in the muscles alleviates pain. Standard size and can be personalized with velcro.

Indications:

- Lower calf muscle stresses
- Protection during sports activities
- Where the lower calf is painful
- The komperasyo

STANDARD CROSS ANKLE SUPPORT



Reference Number:

PF - 917

Sizes:

Universal size

Explanation:

It is made of cotton, breathable soft material and can be used easily in shoes. Provides flexible support to the ankle, relieves edema. Standard size and can be personalized with velcro.

Indications:

- Ankle support
- Gypsum post protection
- Ankle pain

STANDARD SIMPLE ANKLE SUPPORT



STANDARD MALLEOL ANKLE SUPPORT



Reference Number:

PF - 918

Sizes:

Universal size

Explanation:

It is made of neoprene material and it is designed anatomically. Supports the ankle, reducing pain. Neoprene material accelerates blood circulation while maintaining body temperature. Standard size and can be personalized with velcro.

Indications:

- Protection against minor fractures
- 1st degree sprains and strains
- Post-buckling stabilization and protection
- Compression edema after injury

Reference Number:

PF - 919

Sizes:

Universal size

Explanation:

Made of neoprene material. There are pads on both sides of the ankle supporting internal and external malleolus. The eight-fold elastic band provides extra protection. Supports ankle, relieves edema and relieves pain. Standard size and can be personalized with velcro.

Indications:

- Postoperative protection
- After minor fractures
- Protection against sprains
- 1st degree sprains and strains

STANDARD LIGAMENT ANKLE SUPPORT



Reference Number:

PF - 920

Sizes:

Universal size

Explanation:

Made of neoprene material. The inner and outer plastic plates, which provide stabilization around the malleols, are manufactured separately as internal and external. It can be opened completely provides comfort in use. The elastic cross band that wraps in eight form is an extra protection bandage. Supports the ankle with the compression effect, reduces edema and pain. Standard size and can be personalized with velcro.

Indications:

- Protection or support after surgery and injury
- Bond instability and tears
- Protective and edema preventive



SOFT CERVICAL COLLAR



Reference Number:

PF - 570

Sizes:

XS : 24 - 28 cm

Explanation:

It is a soft sponge neck which is designed according to anatomical neck shape and provides high comfort neck support. Relieves headaches by restricting head and neck movements. The removable outer cover is washable.

Indications:

- Cervical spine operations
- Degenerative changes after trauma
- Muscle injuries

Reference Number:

PF - 571

Sizes:

XS : 20 - 25 cm

S : 26 - 31 cm

Explanation:

The inner surface is made of sponge fabric with a towel. Soft and non-sweating material provides comfort in use with velcro adjustment of the shoulders with the desired tension back and provides smooth posture.

Indications:

- Clavicle and sterno-clavicle dislocations and splits.
- Clavicle fractures

ARM SLING NET



ARM SLING SPONGE



Reference Number:

PF - 572

Sizes:

XXS : 23 - 26 cm

XS : 27 - 30 cm

Explanation:

It is made of air perforated fabric. Humerus relaxes the shoulder joint and elbow. Provides shoulder immobilization by fixing the arm to the body. Quality fabric structure provides the comfort of the patient during use.

Indications:

- Post-up immobilization
- Shoulder dislocation and half dislocations
- Sprains
- To prevent the formation of edema in the fracture distal

Reference Number:

PF - 573

Sizes:

XXS : 23 - 26 cm

XS : 27 - 30 cm

Explanation:

The inner surface is made of sponge fabric with a towel. Humerus relaxes the shoulder joint and elbow. Provides shoulder immobilization by fixing the arm to the body. Quality fabric structure provides the comfort of the patient during use.

Indications:

- Post-up immobilization
- Shoulder dislocation and half dislocations
- Sprains
- To prevent the formation of edema in the fracture distal

ELBOW SUPPORT WITH PAD



POSTUREX CORSET



Reference Number:

PF - 574

Sizes:

XS : 24 - 28 cm

Explanation:

Made of neoprene material. An anatomically designed silicone pad is available for compression on muscle groups. Local compression reduces tendon tension and pain.

Indications:

- Painful conditions of the elbow
- Increased joint fluid
- Sprains

Reference Number:

PF - 575

Sizes:

XXS : 55 - 65 cm

XS : 65 - 75 cm

Explanation:

It is made of flexible, breathable, non-sweat-woven fabric. On both sides of the spine there are plastic blades that can be removed and removed. The flexible shoulder straps are crossed under the armpit and crossed to the rear and fastened with velcro straps. Cross bands provide sensory feedback that prevents bad posture disorder by providing postural support. It can be worn easily under clothing.

Indications:

- Posture disorders
- Back and shoulder muscles weakness

CHILD INGUINAL TRUSS (SINGLE SIDED)



Reference Number:

PF - 576

Sizes:

Universal size

Explanation:

Made of flexible cotton fabric. Specially designed anatomic pad and bandage to push back the herniated tissue does not interfere with normal daily activities. It is used to prevent recurrences after surgery and to support healing. It provides the perfect solution for support and comfort to prevent pain and discomfort. The buckled system provides comfort in use.

Indications:

- Single-sided inguinal hernia

CHILD INGUINAL TRUSS (DOUBLE SIDED)



Reference Number:

PF - 577

Sizes:

Universal size

Explanation:

Made of flexible cotton fabric. Specially designed anatomic pad and bandage to push back the herniated tissue does not interfere with normal daily activities. It is used to prevent recurrences after surgery and to support healing. It provides the perfect solution for support and comfort to prevent pain and discomfort. The buckled system provides comfort in use.

Indications:

- Double-sided inguinal hernia

MALLEOL ANKLE SUPPORT



Reference Number:

PF - 578

Sizes:

| | |
|---|------------|
| S | 15 - 17 cm |
| M | 17 - 19 cm |
| L | 19 - 21 cm |

Explanation:

Made of neoprene material. There are pads on both sides of the ankle supporting internal and external malleolus. The eight-fold elastic band provides extra protection. Supports ankle, relieves edema and relieves pain.

Indications:

- Protective and supportive purpose after surgery and injury
- After Gypsum
- In cases of edema, hematoma and effusion
- Severe ligament injuries and instability

SIMPLE ANKLE SUPPORT



Reference Number:

PF - 579

Sizes:

| | |
|---|------------|
| S | 15 - 17 cm |
| M | 17 - 19 cm |
| L | 19 - 21 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. Supports the ankle, reducing pain. Neoprene material accelerates blood circulation while maintaining body temperature.

Indications:

- Protection against minor fractures
- 1st degree buckling and stresses
- Post-buckling stabilization and protection
- Compression in post-injury edema

WRIST SPLINT



Reference Number:

PF - 580

Sizes:

| | |
|---|--------------|
| S | 10,5 - 12 cm |
| M | 12 - 13,5 cm |
| L | 13,5 - 15 cm |

Explanation:

Wrist splint made of neoprene material. It provides immobilization of the wrist in neutral position. It has a bandage that enhances stabilization around the wrist. Neoprene material accelerates blood circulation while maintaining body temperature.

Indications:

- Carpal Tunnel Syndrome
- Immobilization in painful situations of the wrist
- Plaster after fracture
- Injury and after surgery

PATELLA KNEE SUPPORT



Reference Number:

PF - 581

Sizes:

| | |
|---|--------------|
| S | 23 - 25,5 cm |
| M | 23,5 - 28 cm |
| L | 28 - 31 cm |

Explanation:

It is made of neoprene material and it is designed anatomically. The circumference of the patella is open and supported by the pad. It is ideal for patients with pain in the patella and having problems after surgery.

Indications:

- Post-operative protection and rehabilitation
- Protection during sports activities
- Increased joint fluid

DENIS BROWN SPLINT

DENIS BROWN SPLINT (ACTIVE)



Reference Number:

PF - 582

Sizes:

| | |
|------|------------|
| Kısa | 29 - 41 cm |
| Uzun | 42 - 68 cm |

Explanation:

It is the hip abduction orthosis made of aluminum. A pair of boots specially designed for the orthosis is mounted between the two feet. The channel system in which length adjustment can be made is produced. Hip dislocation is used in pesceinovarus deformity.

Indications:

- Hip dislocation

Reference Number:

PF - 588

Sizes:

| | |
|------|------------|
| Kısa | 29 - 41 cm |
| Uzun | 42 - 68 cm |

Explanation:

Splint: Unlike PF-582, it is a dynamic splint, but can perform kicking or crawling movements where the legs are independent of each other.

Indications:

- Hip dislocation

DENIS BROWN BOOT

FREJKA ORTHOSIS



Reference Number:

PF - 583

Sizes:

Number 12
Up to 23 Numbers

Explanation:

Bottom with leather sole, upper section leather with a pull-out, lace-up, front fingertip and heel open to the Denis Browne ateline mountable bottom.

Indications:

Reference Number:

PF - 584

Sizes:

S (1) : 27 - 32 cm

M (2) : 32 - 37 cm

L (3) : 37 - 42 cm

XL (4) : 42 - 47 cm

XXL(5) : 47 - 52 cm

Explanation:

The anatomically designed 12 mm perforated palastosis is made of foam material. Non-hypoallergenic hypoallergenic foam material does not cause skin irritation. X-ray, CT and MRI are compatible. It is used to position the hip joint in flexion and abduction.

Indications:

- Hip dislocation

PAVLIC BANDAGE



Reference Number:

PF - 585

Sizes:

| | |
|---|------------|
| S | 33 - 41 cm |
| M | 41 - 51 cm |
| L | 51 - 61 cm |

Explanation:

Made of polyamide material. Bandage is a dynamic splint that keeps the hip joint in the desired flexion, abduction position and provides a safe range of motion. The tightening of the anterior arch is done by flexing the hip and by tightening the rear belt. It also prevents adduction and extension of the hip. It supports the normal development of the hip joint by allowing the femur to sit in the acetabulum.

Indications:

- Hip dislocation

D.K.C ORTHOSIS



Reference Number:

PF - 586

Sizes:

| | |
|--------|------------|
| S (1) | 27 - 32 cm |
| M (2) | 32 - 37 cm |
| L (3) | 37 - 42 cm |
| XL (4) | 42 - 47 cm |
| XXL(5) | 47 - 52 cm |

Explanation:

Manufactured from anatomically designed thermoplastic material. The inner surface of the towel fabric does not cause skin irritation. It is used to position the hip joint in flexion and abduction.

Indications:

- Hip dislocation

IMMOBILIZER



Referans No:

PF - 587

Beden:

DG-25 : 25 cm

DG-30 : 30 cm

DG-35 : 35 cm

DG-40 : 40 cm

Reference Number:

PF - 587

Sizes:

DG-25 : 25 cm

DG-30 : 30 cm

DG-35 : 35 cm

DG-40 : 40 cm

Açıklamalar:

İç yüzeyi havlu olan süngerli kumaştan üretilmiştir. Dizin ekstansiyonda immobilizasyonunu sağlar. Çıkarılabilir özellikte posterior, medial ve lateralde stabilizasyon sağlayan alüminyum destekleri mevcuttur. Diz yaralanmaları ve ameliyat sonrası için destekleyicidir. Standart beden olup 45 cm, 55 cm, 65 cm uzunlukta boyları mevcuttur. Üç parçalı ped tasarımı sayesinde, kolayca bacak çevresi için ayarlanabilir.

Explanation:

The inner surface is made of sponge fabric with a towel. Provides immobilization of the index extension. Removable aluminum supports are available that provide posterior, medial and lateral stabilization. It is supportive for knee injuries and postoperative surgery. Standard size is 45 cm, 55 cm, 65 cm in length are available. Thanks to the three-piece pad design, it can be easily adjusted for the leg circumference.

Endikasyonlar:

- Menisküs yaralanmalarında ameliyat sonrası
- Akut diz burkulma ve yaralanmaları
- Ameliyat sonrası istirahat, koruma ve immobilizasyona yardımcı
- Ligament yaralanmaları sonrası koruma
- Osteokondral kırıkları
- Diz eklemi travmaları

Indications:

- Meniscal injury after surgery
- Acute knee sprains and injuries
- Helps post-operative rest, protection and immobilization
- Protection after ligament injuries
- Osteochondral fractures
- Knee joint traumas



SILICONE INSOLE

SILICONE INSOLE WITH ARCH SUPPORT



Reference Number:

PF - 444

Sizes:

| | |
|-----|------------|
| XS | 35 - 36 No |
| S | 37 - 38 No |
| M | 39 - 40 No |
| L | 41 - 42 No |
| XL | 43 - 44 No |
| XXL | 45 - 46 No |

Explanation:

It is made of anti-allergic silicone material that does not harm the skin. Specially designed according to anatomical foot shape. The metatarsal heads and the calcaneal epin region, where the body burden is more heavily softened, are specially softened. Reduces pain in feet and calf. It provides balanced standing distribution of body load and reduces stress and stress on joints. Absorbs ground reaction force.

Indications:

- Pes planus
- Plantar fasciitis
- Heel pains
- Metatarsal pain

Reference Number:

PF - 445

Sizes:

| | |
|-----|------------|
| XS | 35 - 36 No |
| S | 37 - 38 No |
| M | 39 - 40 No |
| L | 41 - 42 No |
| XL | 43 - 44 No |
| XXL | 45 - 46 No |

Explanation:

It is made of anti-allergic silicone material that does not harm the skin. Specially designed according to anatomical foot shape. The metatarsal heads and the calcaneal epin region, where the body burden is more heavily softened, are specially softened. Reduces pain in feet and calf. It provides balanced standing distribution of body load and reduces stress and stress on joints. Absorbs ground reaction force.

Indications:

- Pes planus
- Plantar fasciitis
- Heel pains
- Metatarsal pain

SILICONE HEEL CUP



Reference Number:

PF - 446

Sizes:

| | |
|---|------------|
| S | 34 - 37 No |
| M | 38 - 41 No |
| L | 42 - 45 No |

Explanation:

It is made of anti-allergic silicone material that does not harm the skin. Specially designed according to anatomical foot shape. The calcaneal epin region, where the body load is more charged, is specially softened. It relieves the load on the heel and causing pain. Reduces pain in feet and calf.

Indications:

- Heel spurs
- Plantar fasciitis
- Avascular necrosis of calcaneus (sever disease)
- Heel painful conditions

SILICONE TOE SPREADER



Reference Number:

PF - 447

Sizes:

| | |
|---|------------|
| S | 34 - 37 No |
| M | 38 - 41 No |
| L | 42 - 45 No |

Explanation:

It is made of anti-allergic silicone material that does not harm the skin. Suitable for daytime use. Hallux valgus is used to prevent lateral deviation of the thumb.

SILICONE BUNION SHIELD AND TOE SPREADER



Reference Number:

PF - 448

Sizes:

Universal Size

Explanation:

It is made of anti-allergic silicone material that does not harm the skin. Suitable for daytime use. Hallux valgus is used to prevent lateral deviation of the thumb. The anatomic silicone sheath that protects against the friction-preventing trauma on the bunion is used to relieve the pain and pain on the bunion.

CAST SHOE



Reference Number:

PF - 451

Sizes:

| | |
|---|------------|
| S | 33 - 35 No |
| M | 36 - 40 No |
| L | 41 - 45 No |

Explanation:

Top surface of a rigid material which prevents slipping while the bottom part is made of lightweight EVA material. Velcro long waits allow the product to be adjusted to the desired width. In cases where the foot of the patient is cast, it is used to perform normal shoe tasks and facilitates the patient's walk with plaster. It prevents plaster from mechanical damage such as wear and

tear and reduces the impact of the foot.

HALLUX VALGUS NIGHT SPLINT

HALLUX VALGUS SHOE



Reference Number:

PF - 453

Sizes:

Universal size

*There is right and left option.

Explanation:

The inner surface of which is made of sponge cloth towels. By applying a slight pulling force, it prevents lateral deviation of the thumb and progression of the deformity. Suitable for night use.

Indications:

- Hallux valgus deformity

Reference Number:

PF - 452

Sizes:

| | |
|----|------------|
| S | 35 - 37 No |
| M | 38 - 39 No |
| L | 40 - 41 No |
| XL | 42 - 44 No |

Explanation:

Post-op helps early mobilization. It is effective in decreasing the edema in the foot and accelerating the healing process. Long velcro fasteners allow the product to be fitted with the desired width and stiffness. The slope at the bottom is effective in reducing the load on the forefoot.

Indications:

- Foot surgery and trauma

STARFLEX INSOLE WITH ARCH SUPPORT

CORK INSOLE



Reference Number:

PF - 449

Sizes:

Between 35 - 44 Numbers

Explanation:

It is made of special thermo-plastic material. Provides a balanced distribution of the body load, reducing the load on the joints.

Indications:

- Pes planus-planovalgus
- Pes kavus
- Transverse arch needs to be supported

Reference Number:

PF - 450

Sizes:

Between 35 - 44 Numbers

Explanation:

It is made of special cork material. Provides standing balanced distribution of body load, lighten the load on the joints.

Indications:

- Pes planus-planovalgus
- Pes kavus
- Transverse arch needs to be supported



SILICON BREAST PROSTHESIS

BREAST BRA



Reference Number:

PF - 333

Sizes:

15 sizes (60 - 130)

Explanation:

Made of high quality silicone. It is used to provide a natural appearance in cases where breast tissue is taken after mastectomy, to reduce the pain caused by muscle imbalances, and to protect the surgical wound from external impacts and crushes.

Reference Number:

PF - 334

Sizes:

15 sizes (60 - 130)

Explanation:

Prosthetic bra, which is very comfortable to use and stylish in appearance, is specially designed for breast prosthesis. Supports the prosthesis and weight and provides comfortable use. The inner surface is made of cotton fabric that does not irritate the skin.

The shoulder straps are softened for comfortable use. The pocket on both sides allows the prosthesis to be used as right and left.

SILICON EXCERCISE BALLS



SAND WEIGHT



EXCERCISE BAND



Reference Number:

PF - 335

Sizes:

- Red - Hard
- Blue - Medium
- Green - Soft

Explanation:

It is made of ergonomically shaped, high quality, anti-allergic silicone material. Does not leave residue and odor. It takes its former form after use and does not change in the degree of hardness. It provides strengthening of hand, finger and arm muscles, increasing hand skill and mobility, decreasing tension and stress. According to the resistance has 3 different

Reference Number:

PF - 336

Explanation:

Made of neoprene fabric. 1-2-3-4-5 kg options are available. It is used in rehabilitation and strengthening muscle groups in athletes. Suitable for running, rehabilitation or exercise. Adjustable weight bags provide the necessary resistance. It increases the speed of the body's power and metabolism.

Reference Number:

PF - 342

Explanation:

It is very versatile on progressive exercise bands. It allows young and old to do various exercises in the hands of everyone and using it. Suitable for private or general muscle training. There are 4 different color options.




VISCO WAIST SUPPORT PILLOW
COMFORT REF: PF 227



VISCO NECK PILLOW
REF: PF 230



VISCO CLASSIC NECK PILLOW
COMFORT REF: PF 225



VISCO NECK PILLOW COMFORT
REF: PF 223



VISCO TRAVELING PILLOW
REF: PF 231




VISCO ANATOMIC WAIST SUPPORT
PILLOW COMFORT REF: PF 226



VISCO NECK PILLOW
REF: PF 231




VISCO CLASSIC NECK PILLOW
GEL COMFORT REF: PF 224



- Pillows produced with visco elastic memory form technology are sensitive to body temperature, pressure and ambient temperature.
- Products fully support the neck, waist space.
 - Keeps your spine in a correct line during sleep.
 - it allows you to sleep comfortably throughout the night.
 - Reduces your unwanted movements by 80%.
 - Minimizes pressure, regulates blood circulation and removes joint pain.
 - Helps to eliminate problems such as neck stiffness, low back pain and arm numbness.
 - It is hygienic, non-allergic, does not contain carcinogens.

VISKO REFLUX PILLOW
REF: PF 237



SEAT CUSHION (OPEN)
REF: PF 234



WAIST SUPPORT CUSHION
REF: PF 232



VISCO ANATOMIC SITTING CUSHION
COMFORT REF: PF 228



VISCO SEAT CUSHION
REF: PF 236



SEAT CUSHION
REF: PF 235



VISCO SEAT CUSHION COMFORT
REF: PF 229



VISCO WAIST SUPPORT PILLOW
REF: PF 233





WHAT IS VARIS?

It is the main drainage system in the deep veins of the lower extremity, the venous system in the lower extremity. Venous insufficiency of veins is called venous insufficiency, and the veins become apparent.

Although varicose veins usually present as a problem for aesthetics and appearance, they can progress to serious health problems when they are not treated. A method of treating painless, painless varicose veins is the use of varicose veins.

MEDICAL VARISE SOCKS

Puffix Varis socks with world-class pressure profile and knitting quality are a comfortable solution for the need for varicose stockings with its high quality textile surface and robust weaving.

Uygun uygun TS ENV 12718 Varicose Veins are designed according to 127 'standard. CCL1 (18-21 mmHg), CCL2 (23-32 mmHg), CCL3 (34-46 mmHg), pressure and open-closed finger models are available.

The raw materials used are Oeko-Tex 100 certified and anti-allergic. Varicose veins apply a decreasing compression to the calf and upper leg, with the external pressure being highest in the ankle. Thus, they accelerate blood flow by providing a uniform pressure distribution on the leg and on the foot.

INDICATIONS

Low Pressure (18-21 mmHg) Varicose Veins Stockings Medium Pressure (23-32 mmHg)

Payments in certain varicose veins after treatment

Prevention of venous diseases during the recovery of simple ulcers

Post-traumatic swellings leading to expansion and loss of elasticity in the vein

Diabetic and rheumatism type diseases.

Slight varicose veins not accompanied by edema and

early onset of pregnancy,


High Pressure (34-46 mmHg) Varicose Veins

Significant varicose veins accompanied by edema

In the recovery period of severe and recurrent ulcers (after the ulcer has healed)







In severe structural and post-thrombotic venous insufficiencies

Lymph edema treatment



Model

CCL 1 CCL 2 CCL 3
18 - 21 mmHg 23 - 32 mmHg 34 - 46 mmHg

| | | | | |
|---|---|--------|--------|--------|
|  | A-D Diz Altı Kapalı Parmak | PV-101 | PV-107 | PV-113 |
|  | A-D Diz Altı Açık Parmak | PV-100 | PV-106 | PV-112 |
|  | A-G Diz Üstü Kapalı Parmak | PV-103 | PV-109 | PV-115 |
|  | A-G Diz Üstü Açık Parmak | PV-102 | PV-108 | PV-114 |
|  | A-T Kilotlu Kapalı Parmak | PV-105 | PV-111 | PV-117 |
|  | A-T Kilotlu Açık Parmak | PV-104 | PV-110 | PV-116 |

Over The Knee Varis Sock

| Beden | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
|--------------------------|----------|---------|---------|---------|---------|---------|---------|
| Ayak Bileđi Çevresi (cb) | 18-20 cm | 20-22cm | 22-24cm | 24-26cm | 26-28cm | 28-30cm | 30-32cm |

Under The Knees Varis Sock

| Beden | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
|--------------------------|----------|---------|---------|---------|---------|---------|---------|
| Ayak Bileđi Çevresi (cb) | 18-20 cm | 20-22cm | 22-24cm | 24-26cm | 26-28cm | 28-30cm | 30-32cm |
| Alt Baldır Çevresi (cC) | 18-33 cm | 30-37cm | 33-40cm | 35-43cm | 37-45cm | 39-46cm | 41-48cm |
| Üst Baldır Çevresi (cF) | 37-47 cm | 38-52cm | 43-58cm | 46-61cm | 50-66cm | 53-69cm | 56-71cm |

Pantyhose Stockings

| Beden | 1 | 2 | 3 | 4 | 5 | 6 | 7 |
|--------------------------|----------|---------|---------|---------|---------|---------|---------|
| Ayak Bileđi Çevresi (cb) | 18-20 cm | 20-22cm | 22-24cm | 24-26cm | 26-28cm | 28-30cm | 30-32cm |
| Alt Baldır Çevresi (cC) | 18-33 cm | 30-37cm | 33-40cm | 35-43cm | 37-45cm | 39-46cm | 41-48cm |
| Kasık Çevresi (cG) | 43-57 cm | 46-60cm | 49-67cm | 52-72cm | 55-76cm | 60-79cm | 64-80cm |

ANTI-EMBOLI SOCKS

Anti-embolism socks accelerate blood flow and reduce coagulation by applying constant pressure to the legs. The external pressure exerts a decreasing compression towards the calf and upper leg, the highest in the ankle. Puffix Anti-Embolism Stockings si oz TS ENV 12719 Medical Thrombosis Prophylaxis Socks Anti Anti is designed in accordance with the standard. There are 3 models, including knees, laptops and notebooks. There are anti-slip silicone bands on top of the sock. Thanks to its special latex-free weaving, it has a high elasticity.

INDICATIONS



Increases blood flow rate by providing graded printing.
It reduces the risk of deep vein thrombosis in immobilized patients before and after surgery.
Increases blood flow in other vein disorders such as chronic venous stasis.

CONTRAINDICATIONS

Pulmonary edema due to severe arteriosclerosis or other ischemic vascular disease
Cellulitis Skin inflammation Gangrene Fascitis, Recently skin graffiti, Paniculitis
Pulmonary edema due to severe edema or subjective heart disease in the legs
Extreme deformity of the leg

Model

CCL 13-18mmHg

| | | | |
|--|-----------------|--------|---------|
|  | A-D Diz Altı | PV-113 | |
|  | A-G Diz Üstü | PV-114 | PV-115K |

Under The Knees - Anti Emboli Socks

| Beden | XS | S | M | L | XL |
|---|----------------|-----------------|------------------|--------------------|-----------------|
| Ayak Bileği Çevresi (cB) Topuk Rengi | 18-20cm Bej | 20-23cm Sarı | 23-26cm Yeşil | 26-29cm Kırmızı | 29-32cm Mavi |

Over The Knees - Anti Emboli Socks

| Beden | XS | S | M | L | XL |
|---|--------------------------------------|---------------------------------------|--|--|---------------------------------------|
| Ayak Bileği Çevresi (cB) Alt Baldır Çevresi (cC) Üst Baldır Çevresi (cF) Topuk Rengi | 18-20cm 27-34cm 27-34cm Bej | 20-23cm 31-37cm 31-37cm Sarı | 23-26cm 35-40cm 35-40cm Yeşil | 26-29cm 38-44cm 38-44cm Kırmızı | 29-32cm 41-48cm 41-48cm Mavi |

puffix
Yetiřkin Hasta Bezi
Adult Diaper



puffix
Yatak Koruyucu Örtü
Bed Sheet



Markalarıdır



YETKİLİ BAYİ

Tüm görsellerin telif hakkı Atılım Sağlık A.Ş. ne ait olup izinsiz kullanılamaz.
Katalogdaki ürünlerin model, renk ve ebatlarında değişiklik olabilir.

Adres: Merkez Mh. Değerimenbahçe Cd. Kavak Sk. No:9 Ağaçlı Silo Binası K:1 Yenibosna - Bahçelievler - İstanbul
Tel: 0212 667 77 87 **Faks:** 0212 667 77 07 www.atilimsaglik.net